

# *Big Lifter*<sup>TM</sup>

**BRUNO**  
for your independence



*Superior Performance.*

*Amazingly Compact.*

*Easy To Use.*

**VSL-570 Vehicle Lift**

## *Versatile Manual Rotate Lift*

The VSL-570 lifts up to 400 lb/  
181 kg at the touch of a button.  
Graceful, automotive styling and  
simple operation make this product  
a great solution to transport a scooter  
or powerchair in your vehicle.



**Fold-Down Lift Head (requires optional quick release pins)**

### *Value Added Features...*

- Hundreds of vehicle specific mounting bases provide clean and precise installation into minivans, SUVs, CUVs and full-size vans
- Fold-Down Lift Head provides greater visibility and storage space. Optional quick release pins make this a snap.
- Upgradeable to power rotate function
- Lift design allows flexibility to raise and lower your mobility device behind, at the corner, or at the side of your vehicle

# *Innovation... for your Independence & Convenience.*

## **KEY SPECIFICATIONS**

	<b>VSL-570</b>
Maximum Weight Capacity	400 lb/181 kg
Hand-Held Pendant Controls	yes
Adjustable Height	yes
Adjustable Length/Angle Lift Head	yes
Manual Rotation	yes
Upgradeable To Power Rotate	yes
Power Raise/Lower	yes
Fold-Down Lift Head Option	yes
Predominately Stainless Steel Hardware	yes
Eligible for Auto Manufacturer Rebate	yes
Exterior Option for Pickup Trucks	yes
3 Year Warranty	yes

**IMPORTANT:** Contact your factory trained Bruno dealer for specific details on your vehicle!



**1.** Connect docking device...



**2.** Raise mobility device by pressing button...



**3.** Manually rotate mobility device into vehicle.

**Represented by:**

**BRUNO INDEPENDENT LIVING AIDS, INC.®**

1780 Executive Drive, P.O. Box 84, Oconomowoc, WI 53066  
www.bruno.com

