



# Veigel right hand controls



## The Classic Hand Control 4100CL

Our Classic hand controls are made for simple and reliable operation. To accelerate, simply turn the ergonomic handle clockwise. To brake push the unit slightly forward. The brake can be locked while waiting in traffic. The ergonomic handle can be adjusted to accommodate individual driver needs. The unit can be folded out of the way for ambulatory drivers.



## Commander for Classic Hand Control 4100CLC

Our Classic hand control can be upgraded to add a five button switch to control secondary functions like the turn signals, wipers, horn, etc... The commander upgrade comes complete with a relay box, wiring diagram and wiring harness.



## The Compact Hand Control 4200CP

Proven technology, sleek and elegant design - that is how the Compact hand control is made. The function is simple: to accelerate, pull the ergonomically shaped handle backwards, a slight push in the forward direction operates the brake. The brake can be locked at the push of a button, while waiting in traffic. The handle is shaped in such a way that the hand can rest on it while driving. This means relaxed and effortless driving, even on long distances. The hand control can be folded out of the way when not in use.



## Commander for Compact Hand Control 4200CPC

Our Compact hand control can be upgraded to add a five button switch to control secondary functions like the turn signals, wipers, horn, etc... The commander upgrade comes complete with a relay box, wiring diagram and wiring harness.



## The Basic Hand Control 4300B

Our Basic hand controls are simplified, no frills version of the Compact hand control. Offers first rate quality at a competitive price.