



# Fenton Mobility

PRODUCTS

## 2008 Ford Econoline Wagon XLT

Vin	1FMNE11L58DB50689
Miles	34998
Stock	B50689
Engine	Other 5.4L EFI V8 ENGINE
Warranty	
E. Color	
I. Color	
Status	sold
Trans.	Automatic
Drive	
Body	Minivan/Van
<b>Price:</b>	<b>\$32,995</b>



Price After All Rebates & Discounts. Does Not Include Tax, License and/or Other Fees. Please Call For Exact Pricing. Pricing Subject To Change. The information on this page is subject to change. Please call in advance of coming out to confirm availability, options & equipment. The list of data on this page does not necessarily reflect what is actually on the vehicle. **Always confirm options before purchasing.**

### Options:

#### Ability Van Options:

- ✓ 60 Inch Doorway Opening
- ✓ Ability Series Door Struts
- ✓ Auxiliary Stainless Steel Cables And Loops
- ✓ Ford Chassis Equipment Super Cargo Vans
- ✓ Full Length Body Graphics Upper And Lower
- ✓ NHTSA Compliant Interlock
- ✓ Rear Air Conditioning
- ✓ Seamless Headliner With Recessed Dome Lights
- ✓ Ability 20 Inch Aerodynamic Roof With FMVSS Certified
- ✓ Ability Series Side Walls Providing Maximum Space
- ✓ Back Up Alarm
- ✓ Front A/C Heat
- ✓ Heat Passenger Side Mount With Rubber Fluid Lines
- ✓ NHTSA Compliant Lift Install Kit
- ✓ Recessed L Series Tie-down Track
- ✓ Ability Series Door Aluminum Door Grabbers
- ✓ Ability Series Stainless Steel Raised Doors
- ✓ Colored Body
- ✓ Full Body, Privacy Glass
- ✓ Medium Flint Vinyl Front Bucket Seat
- ✓ Overhead Storage With T Latch
- ✓ Seamless Flooring

### Description:

ABILITY PT The Ability PT conversion offers all the accessibility without all the added cost. We start with a Ford Passenger van and add Stainless Steel raised side doors giving consumers 56 of clear passage through the side door. Then we add a sleek raised roof with a certified roll cage offering plenty of interior height. Also included is a remounted rear factory bench seat over head air conditioning center section wheelchair flooring and tie-down tracks. A Braun Millennium lift is added finished off with automatic door struts which automatically keep the doors out of your way during lift operation. You wont find a better price for this vehicle anywhere.

### Vehicle Location

Contact Fenton Mobility Products Inc.  
Phone: 716-633-7016  
Address:  
2851 Broadway Suite 900  
Buffalo, NY 14227  
<http://www.fentonmobility.com>

### About Location

From the city take 33 West just past the Airport on your right. From the west take Transit road to Genesee Street. We are on the left just before the Airport. Directly accross from the Noco/Holiday Inn

### Sales Hours:

Monday	8:30am	To: 5:00pm
Tuesday	8:30am	To: 5:00pm
Wednesday	8:30am	To: 5:00pm
Thursday	8:30am	To: 5:00pm
Friday	8:30am	To: 5:00pm
Saturday	By	To: Appt.
Sunday	Closed	To: Closed

### About Us:

Fenton Mobility is an industry leader in the modification and adaptation of wheelchair vans & handicap vans. Our trained and professional staff is available to assist you with your adaptive equipment needs, from sales to service, consumer handicap vans to commercial wheelchair vans. Our customers are provided with an extensive personal consultation, fittings and overview of all the necessary mobility equipment and van modifications that will allow freedom and mobility for years to come. Special Mobility Financing is available. [Read More Here](#)